



Pricing for *TUCBE* instrument and Accessories.

Standard price list March 2008

<p>TUCBE includes, NeCu laser at 248.6nm, Basic 6 channel Fluorescence detector, User defined filter set, Standard(280, 300, 320, 340, 360, 380nm all 20nm FWHM), PMT response 180-650nm, pulse by pulse input energy meter, LabView control, analysis and system calibration software, operation and application manuals, JAUS Compliant 1 year warranty on complete system.</p>		\$ 40,000
<p>Photo Multiplier Tube option. (pmt link)</p>	NEW	\$700
<p>PMT spectral response 185nm- 850nm -4</p>	Replace	\$160
<p>PMT spectral response 185nm- 650nm (included) -6</p>		
<p>Computer, It is recommend that you purchase the laptop computer loaded with XP, LabView exe, all the drivers for the laser and monochrometer, Excel, All tested and with original manufacturer warrantee. Usually a Dell, call for specs.</p>	*Price may vary depending on spec.	*\$1500
<p>Pelican hard case with custom foam, TUCBE Peripherals. Roll-away for easy transportation.</p>		\$400
<p>Battery operating option: Use the system in a remote location. Provide is a Battery system that can be used in conjunction with the TUCBE system giving it a dynamic operating range of >4 hours continuous use. Includes two Batteries and 1 Quad charger system.</p>		\$3,100